

Ultrafast Plasma Protein Binding Analysis

Using the Agilent RapidFire high-throughput mass spectrometry system and accurate mass Q-TOF

Authors

Kari E. Schlicht and
Vaughn P. Miller
Agilent Technologies, Inc.

Abstract

There has been increased demand for higher throughput *in vitro* ADME analyses to enable earlier optimization of drug candidates. One analysis that is key to assessing both the pharmacokinetics and pharmacodynamics of a drug candidate is plasma protein binding (PPB). A single, ultrafast method for analysis using the Agilent RapidFire High-Throughput Mass Spectrometry system and an Agilent 6550 iFunnel Q-TOF Mass Spectrometer was developed. More than 40 compounds were subjected to rapid equilibrium dialysis, then analyzed using a 9 second/sample SPE/TOF/MS method as well as by LC/MS/MS. Excellent correlation was determined between the two analytical methodologies with an R^2 value of 0.977. However, the SPE/TOF/MS method provided 20 times greater throughput compared to LC/MS/MS with a capacity of 400 samples/hour.

Introduction

During plasma protein binding, the fraction of drug that is bound to plasma proteins is no longer available to interact with the biological target or site of action. The remaining unbound and bioavailable fraction, therefore, is an important factor in pharmacokinetic properties including distribution and clearance. LC/MS/MS has become the method of choice for determining PPB due to its selectivity and sensitivity. However, disadvantages include the need to develop individual MRM methods for every compound and analysis times of several minutes per sample. In this study, these time-consuming disadvantages were overcome by eliminating the need for individualized compound method development and dramatically increasing the speed of analysis. The ability to develop a single, ultrafast method of analysis using the RapidFire High-Throughput Mass Spectrometry system and a 6550 Q-TOF Mass Spectrometer was investigated.

Experimental

Analytical systems

An Agilent 1260 Infinity Binary LC system was coupled with an Agilent 6460 Triple Quadrupole Mass Spectrometer and used to analyze each sample. Samples were analyzed at a rate of 3 minutes per sample using the conditions shown in Table 1. Individual compound multiple reaction monitoring (MRM) methods were optimized off line by infusion using a syringe pump from KD Scientific, Holliston, MA.

Identical samples were also analyzed by a SPE/TOF/MS system, consisting of a RapidFire 360 High-throughput Mass Spectrometry system, and a 6550 iFunnel Q-TOF Mass Spectrometer. Samples were analyzed at a rate of 9 seconds or less per sample using the conditions shown in Table 2.

Table 1. LC-triple quadrupole instrument conditions.

| LC Conditions | |
|------------------------------|--|
| Analytical Column | Agilent Poroshell 120 EC-C18, 2.1 mm × 50 mm, 2.7 μm (p/n 699775-902) |
| Column Temperature | 25 °C |
| Injection Volume | 5 μL |
| Mobile Phase | A) 0.1 % formic acid in water B) 0.1 % formic acid in acetonitrile |
| Gradient | 10 % B 0.0 minutes 10 % B 0.2 minutes 95 % B 0.8 minutes 95 % B 2 minutes |
| Run Time | 2.0 minutes |
| Flow Rate | 0.6 mL/min |
| Post Time | 1 minute |
| Triple Quadrupole Conditions | |
| Acquisition Mode | Agilent JetStream, positive ionization, MRM |
| Sheath Gas Temperature | 350 °C |
| Sheath Gas Flow Rate | 11 L/min |
| Drying Gas Temperature | 300 °C |
| Drying Gas Flow Rate | 10 L/min |
| Nebulizer | 35 psi |
| Nozzle Voltage | 500 V |
| Capillary Voltage | 3,500 V |

Table 2. RapidFire-Q-TOF instrument conditions.

| RapidFire Conditions | |
|------------------------|--|
| Buffer A | Water with 0.09 % formic acid, 0.01 % trifluoroacetic acid |
| Buffer B | Acetonitrile with 0.09 % formic acid, 0.01 % trifluoroacetic acid |
| Injection Volume | 10 μL |
| SPE Cartridge | Agilent RapidFire cartridge C (reversed-phase C18 chemistry, p/n G9205A) |
| RF State 1 | Sip sensor |
| RF State 2 | 3,000 ms |
| RF State 3 | 3,000 ms |
| RF State 4 | 500 ms |
| MS Conditions | |
| Ion Mode | Dual AJS ESI, positive ion polarity |
| Drying Gas Temperature | 200 °C |
| Drying Gas Flow Rate | 18 L/min |
| Nebulizer | 30 psi |
| Sheath Gas Temperature | 350 °C |
| Sheath Gas Flow Rate | 12 L/min |
| Nozzle Voltage | 300 V |
| Capillary Voltage | 5,000 V |
| Fragmentor Voltage | 175 V |
| Skimmer Voltage | 65 V |
| OCT1 RF Vpp | 750 V |
| RF Voltages | High pressure funnel: 180 V Low pressure funnel: 80 V |
| Acquisition Parameters | 2 GHz dynamic range MS mode 125 to 1,000 m/z acquisition 5 spectra/s |

Chemicals and reagents

Rapid Equilibrium Dialysis (RED) device plates and inserts were purchased from Thermo Scientific (Rockford, IL). Pooled human plasma was purchased from Bioreclamation Inc. (Westbury, NY). Phosphate buffered saline (PBS), test compounds, and solvents were purchased from Sigma-Aldrich (St. Louis, MO).

Sample preparation

Human plasma (100 μ L) was spiked with individual test compounds at a concentration of 5 μ M and inserted into the sample chamber of the RED device. PBS (300 μ L) was inserted into the buffer chamber of the RED device. Triplicate assays were conducted for each test compound. The test plates were sealed and incubated for 5 hours at 37 $^{\circ}$ C on an orbital shaker. After incubation, 50 μ L from each sample chamber was transferred to a well of a 96-well plate containing 50 μ L of PBS. The corresponding 50 μ L from each buffer chamber was transferred to a well of the 96-well plate containing 50 μ L of human plasma. Cold acetonitrile (200 μ L) containing 0.1% formic acid and internal standard (1 μ M bucetin) was added to each sample to precipitate proteins. The samples were centrifuged, and the resulting supernatants were diluted 1:1 with water and transferred to a 96-well analysis plate.

Data analysis

Agilent MassHunter Quantitative Analysis software (B.05.00) was used to integrate samples acquired from LC/MS/MS analysis. Samples acquired from SPE/TOF/MS were extracted by exact mass using a 10 ppm window and integrated using RapidFire Integrator Software. All integrated values were normalized by their corresponding

internal standard values. Percentage bound was calculated for each test compound using the following equation:

$$\% \text{ bound} = 100 - \left[\left(\frac{\text{average buffer}}{\text{average plasma}} \right) \times 100 \right]$$

Results and discussion

More than 40 different test compounds with diverse chemical properties including a xLogP range of -0.4 to 7.1 and molecular weight range of 160 to 733 g/mol were subjected to rapid equilibrium dialysis. Identical samples were analyzed by SPE/TOF/MS and by LC/MS/MS. A small, chemically diverse, subset of the test compounds was used to develop a single generic SPE/TOF/MS method. Sample-to-sample rates for this method were 9 seconds or less, providing a throughput of 400 samples an hour. This rate is 20 times greater

than the throughput of the LC/MS/MS method. In addition, the LC/MS/MS method with gradient chromatography required MRM method development before analysis for every compound. This bottleneck in method development was eliminated by the use of mass accuracy provided from the high resolution of the 6550 Q-TOF.

The analytical results for the SPE/TOF/MS method were found to be comparable to traditional LC/MS/MS despite the differences in methodology. The percentage bound determined for each compound by both analyses is shown in Table 3. There was excellent correlation between the two methods with a R^2 value of 0.977 (Figure 1). There was also good agreement for both methods with previously reported values obtained using RED devices in the literature.¹⁻³

Table 3. Human plasma protein binding values.

| Compound | LC-Triple Quadrupole, % Bound | RapidFire-TOF, % Bound | Difference | RED Lit. |
|------------------|-------------------------------|------------------------|------------|----------|
| Amitriptyline | 98.0 | 97.1 | 0.90 | 98.9 |
| Amodiaquine | 99.0 | 94.9 | 4.15 | |
| Amoxapine | 92.3 | 93.6 | -1.26 | |
| Buspirone | 89.5 | 86.5 | 2.97 | |
| Carbutamide | 91.4 | 91.2 | 0.19 | |
| Chlorpromazine | 99.9 | 99.8 | 0.12 | |
| Cinnarizine | 99.6 | 99.3 | 0.25 | |
| Desipramine | 90.7 | 91.9 | -1.23 | |
| Dextromethorphan | 69.7 | 70.4 | -0.74 | |
| Diclofenac | 99.7 | 99.7 | -0.03 | |
| Diltiazem | 87.9 | 82.6 | 5.29 | |
| Diphenhydramine | 84.8 | 84.0 | 0.79 | |
| Erythromycin | 82.8 | 80.4 | 2.41 | 79.1 |
| Fluconazole | 20.7 | 17.6 | 3.05 | 21.2 |
| Fluphenazine | 99.4 | 98.8 | 0.62 | 95.7 |
| Imipramine | 94.6 | 93.6 | 0.96 | 93.4 |
| Ketoconazole | 99.4 | 99.3 | 0.03 | |
| Lansoprazole | 99.0 | 98.9 | 0.13 | |
| Levofloxacin | 46.1 | 50.5 | -4.31 | 60.3 |
| Metoprolol | 24.3 | 16.7 | 7.60 | 3.5 |
| Nadolol | 23.0 | 17.7 | 5.29 | 25.3 |
| Nicardipine | 99.9 | 99.8 | 0.05 | |
| Nicotine | 54.4 | 42.5 | 11.94 | |

Conclusion

An ultrafast SPE/TOF/MS method was developed to analyze plasma protein binding assay samples at a rate of 9 seconds per sample with analytical results comparable to those determined by LC/MS/MS. In addition to a 400 sample/hour throughput, this method provided increased efficiency by eliminating the need for offline MRM method development. This ultrafast, high-resolution approach may be useful for the analysis of other *in vitro* ADME assays.

References

1. van Liempd, S. *et al.* Development and Validation of a Higher-Throughput Equilibrium Dialysis Assay for Plasma Protein Binding. *Journal of the Association for Laboratory Automation* **2011**, *16*(1), 56–67.
2. Waters, N. J. *et al.* Validation of a Rapid Equilibrium Dialysis Approach for the Measurement of Plasma Protein Binding. *Journal of Pharmaceutical Sciences* **2008**, *97*(10), 4586–95.
3. Shanler, M. S. *et al.* Validation of an Automated High Throughput Plasma Protein Binding Assay. *BD Biosciences application note #474*, **2009**.

| Compound | LC-Triple Quadrupole, % Bound | RapidFire-TOF, % Bound | Difference | RED Lit. |
|--------------|-------------------------------|------------------------|------------|-----------|
| Nifedipine | 97.6 | 91.6 | 5.91 | |
| Nimodipine | 98.8 | 98.4 | 0.37 | |
| Nizatidine | 33.4 | 18.0 | 15.39 | |
| Norfloxacin | 55.4 | 47.0 | 8.36 | |
| Phenacetin | 37.2 | 44.9 | -7.70 | |
| Promazine | 97.7 | 97.9 | -0.24 | |
| Promethazine | 99.1 | 98.6 | 0.47 | 98.1 |
| Propafenone | 95.9 | 94.0 | 1.91 | |
| Propranolol | 90.3 | 87.2 | 3.09 | 80.2-92.6 |
| Pyrilamine | 72.3 | 69.8 | 2.51 | |
| Quinidine | 88.5 | 85.7 | 2.85 | |
| Tacrine | 74.2 | 63.1 | 11.12 | 63.3 |
| Tamoxifen | 100.0 | 99.5 | 0.46 | |
| Terfenadine | 95.7 | 99.9 | -4.20 | 99.5 |
| Testosterone | 95.2 | 97.1 | -1.92 | |
| Thioridazine | 100.0 | 99.8 | 0.19 | 99.9 |
| Ticlopidine | 99.9 | 99.3 | 0.54 | |
| Tolbutamide | 96.3 | 97.7 | -1.37 | 95.6 |
| Tripolidine | 91.3 | 93.1 | -1.75 | |
| Verapamil | 96.3 | 94.2 | 2.08 | 90.6-94.8 |
| Warfarin | 99.4 | 99.5 | -0.08 | 99.5 |




Figure 1. Comparison of human PPB values from LC-triple quadrupole versus RapidFire-TOF.

www.agilent.com

This information is subject to change without notice.

© Agilent Technologies, Inc. 2013, 2022
 Printed in the USA, April 14, 2022
 5991-2497EN