



HIGH CAPACITY LABORATORY GLASSWARE WASHER (365 LITRES)

Professional line - large dimension glassware washers

- washing up to 95°C - antibacterial thermal disinfection
- electronic programming: 20 set programs + 10 free programs
- 3 microprocessors programmable electronic control
- serial port RS232 for connection to external printer or PC
- 128x64 pixel graphics display for different languages messages
- active thermodynamic drying
- peristaltic pump for liquid detergent
- peristaltic pump for liquid neutralizer
- possibility of installing up to 3 additional peristaltic pumps
- washing chamber and internal door made in AISI 316L stainless steel
- all internal piping made in AISI 304 stainless steel
- external panels made in AISI 304 stainless steel
- detergent cabinet made in AISI 304 stainless steel for 4 5l cans
- washing up to 3 levels
- external dimensions LxWxH (mm) 900x800x1900
- internal working dimensions LxWxH (mm) 670x650x835
- 2 independent washing pumps (one for chamber, one for trolleys)
- high efficiency drain pump
- built-in boiler for fast cycles preheating water
- cold water, hot water and demineralized water connections
- electronic door lock system for user safety
- electrical connection 3P+N 400V 50Hz 18kW max
- 60Hz version available
- without trolleys

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350

Phone +39 0522 8211 - Fax +39 0522 821453 - service.instruments@smeg.it